



www.bmtmedikal.com.tr



Ürün Katalođu

“sađlık teknolojisindeki son nokta”




www.bmtmedikal.com.tr

ÜRÜNLERİMİZ

- HASTABAŞI MONİTÖRLERİ
- EKG CİHAZLARI
- PULSEOKSİMETRELER
- FETAL MONİTÖRLER
- ENJEKTÖR POMPALARI
- İNFÜZYON POMPALARI
- KÜVÖZ VE FOTOTERAPİ
- DEFİBRİLATÖRLER
- İNVAZİV KABLolar, İBP KABLolar
- ELEKTRO KOTER AKSESUARLARI
- İSİ PROBLARI
- EKG KABLO , EKG LEAD
- SPO2 PROB
- SPO2 KABLO
- NİBP MANŞON , ADAPTOR ve HORTUMLARI
- OKSİJEN SENSÖRLER, FLOW SENSÖRLER
- VENTİLATÖRLER
- ANESTEZİ CİHAZLARI
- NEMLENDİRİCİLER
- VENTİLATÖR AKSESUARLARI

Perinatal Care Accessories



A	ThermoTrace Core, 5 pcs. Skin temperature probes (yellow) to measure the near-to-core skin temperature, angled.	MX11000
B	ThermoTrace Peripheral, 5 pcs. For peripheral skin temperature measurement (ThermoMonitoring), angled, white.	MX11001
C	ThermoPad, 50 pcs. Adhesive reflecting pads for attachment of the skin temperature probes.	MX11002
D	Bronchial aspirator, complete With ejector up to -0.9 bar, fitted beneath the bed area. Portable secretion jar set, connection to central O ₂ or medical air system. The ejector generates infinitely adjustable vacuum up to -0.9 bar.	2M21189
E	Ventilation hose holder BabyTherm® Ventilation hose support mounted at the Babytherm side walls.	2M21191
F	Deposit plate Table for BabyTherm 8000 WB.	2M21158
G	Handle, complete BT 8000 WB Complete handle for BabyTherm 8000 WB.	2M21161

Not all articles are available worldwide.

Perinatal Care Accessories

A Ventilation hose holder Caleo®, complete 8411075

Hose support for Incubators is affixed directly to the bed base and therefore holds the hoses in the given position even if the bed is raised or lowered or is removed from the incubator, 40 cm (15.8 in) long.

B Air filter Caleo, 20 pcs. MX17015

Fresh air filter for Caleo.

C Infusion holder, complete 2M21514





Infusion support for attachment to 38 mm (1.5 in) pole, with 3 hooks for infusion bags or infusion bottles, e.g. for the option servocontrolled humidity. Necessary accessories:
- 1 x pole 38 mm (1.5 in)/600 (1 x 2M50691).

D Notebook holder 2M22171

Holder on which to place a notebook. Can be attached to standard rail on the wall side.

E O₂ Sensor OxyTrace MX01050

Oxygen probe requiring no calibration for one year, useful life two years.



Perinatal Care

Perinatal Care Accessories



A 2.2 Accessories for Babylog® 8000/8000 plus


2.2.1 Fisher & Paykel hose sets Usage with MR 850.

B Fisher & Paykel hose set 8411041
Hose set with reusable hoses of silicone.


Consists of:

1	Breathing hose, P, 1.1 m (3.6 ft) Silicone breathing hose (pediatric) with 10/22 (0.4/0.9) standard sleeves.	2165651
2	Silicone hose, P, 1 m (3.3 ft)	8403080
3	Silicone hose, P, 35 cm (13.8 in)	8403070
4	Test lung, complete (bellows)	8409742
5	Watertrap 10 Watertrap for neonatal hoses.	8409627
6	F&P Adapter	8409625
7	Double cone 11A	8409897
8	Connectors for tracheal tube, size 11	M19351
9	T-piece adapter Connector with hole for temperature sensor, e.g. for F&P.	8411044
10	Hose heater, 1.1 m (3.6 ft)	8411045
11	Corrugated hose, flexible, 20 cm (7.9 in)	8410709
12	Silicone hose, P, 50 cm (19.7 in)	8418458

Not all articles are available worldwide.



Perinatal Care Accessories



A Fisher & Paykel HFV hose set 8411153

Hose set HFV F&P with autoclavable, lightweight hoses of Hytrel with shorter lifespan and specifically recommended for use of Babylog® 8000 HFV.

Consists of:


1	Hose, lightweight 70 cm (27.6 in)	8418404
2	Hose, lightweight 65 cm (25.6 in)	8410814
3	Test lung, complete (bellows)	8409742
4	Adapter 900 MR 533	8418407
5	Adapter 900 MR 830	8418455
6	Connector for tracheal tube, size 11	M19351
7	Watertrap 10 Watertrap for neonatal hoses.	8409627
8	F&P Adapter	8409625
9	Hose heater	8411154
10	Hose, lightweight 40 cm (15.7 in)	8410816
11	T-piece adapter For F&P temperature sensor.	8418405

Not all articles are available worldwide.

Perinatal Care Accessories

2.5.2 Y-pieces, connectors and cables

A	Y-piece, child, 90° For masks with ISO connectors and standard connectors, with 10 mm (0.4 in) socket connectors and suction opening.	8403075
	Y-piece, child, straight Reusable Y-piece for pediatric hose sets.	M27077
	Y-piece, 90° Child Y-piece.	8403077
B	LuerLock angled probe for neonates LuerLock connector for child Y-piece (M26700).	M34568
C	Neonatal flow sensor Y-piece Y-piece with flow sensor.	8410185
C	Neonatal flow sensor insert, 5 pcs. Insert for flow sensors for neonates, for Babylog 8000 series. 5 pcs.	8410179
D	Flow sensor cable Connecting cable for Dräger flow sensors of Babylog series.	8409626
E	Neonatal flow sensor ISO 15 Straight flow sensor.	8411130
F	Analog cable, 3 pcs. Analog cable for Babylog 8000.	8306487



Not all articles are available worldwide.

TAMİR BAKIM ONARIM

0216 518 65 45-46

01010101000101010001010100011010010011

01010101000101010001010100011010010011

www.bmtmedikal.com.tr

BMT Medikal olarak her marka ve model tıbbi cihazların teknik servis ve bakımı yapılır.



www.bmtmedikal.com.tr



*Bu katalogta bahsedilen ürünler, logolar ve şirket isimleri başka firmaların ticari markaları olabilir. Bu durum yukarıda adı geçen markaların Bmt Medikal ile herhangi bir yakınlığı veya onayı olarak algılanmamalıdır. Bu dokümanda verilen bilgiler dikkatlice kontrol edilmiş ve doğruluğuna inanılmıştır. Sürekli ürün gelişimi adına, Bmt Medikal bu dokümanda ve tarif edilen üründe uyarı ve bilgilendirme yapmadan değişiklikler ve gelişmeler yapma hakkını saklı tutar.



www.bmtmedikal.com.tr



BMT MEDİKAL SİST. VE ÜRÜNL. SAN. ve DIŞ TİC. LTD.ŞTİ.
Aydınevler Mah. İskele Sk. No:13/1 Maltepe-İstanbul - TURKEY
Tel: 0216 518 65 45 - 46 Fax: 0216 518 65 47
e-mail: info@bmtmedical.com web: www.bmtmedical.com